


Vascular Intervention // **Peripheral**
PTA Balloon Catheter/0.035"/OTW

Passeo[®]-35 Xeo

Low profile PTA balloon with excellent deliverability in a wide range of sizes.



Improved crossability



Excellent deliverability



Low profile, wide range of sizes




BIOTRONIK
excellence for life

Passeo-35 Xeo

Low profile PTA balloon
with excellent deliverability
in a wide range of sizes.

Improved crossability

Up to 0.16 mm lower tip entry profile compared to leading competitors.¹





Advanced tip profile

The low tip entry profile and smooth tapered tip of Passeo-35 Xeo facilitates crossing.




Lowest crossing profile

Up to 0.22 mm lower crossing profile compared to leading competitors.¹




Low entry profile
and smooth taper
for better crossing





Excellent deliverability

The Passeo-35 Xeo catheter offers excellent trackability and pushability compared to competitive 0.035" PTA balloon catheters.²




Low profile, wide range of sizes

Passeo-35 Xeo features a low profile design with a full 5F compatibility up to 7.0 mm balloon diameter and 250 mm balloon length, 6F up to 10.0 mm and 7F up to 12.0 mm diameter.



Passeo-35 Xeo

Low profile PTA balloon designed to cross



Passeo-35 Xeo

The Passeo-35 Xeo is indicated to dilate stenosis in the iliac*, femoral, popliteal and infrapopliteal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. Passeo-35 Xeo is also recommended for post-dilation of balloon expandable and self-expanding stents in the peripheral vasculature.**

Vascular
Intervention
Peripheral



Technical Data - Balloon Catheter

Catheter type	OTW	Balloon markers	2 swaged markers
Recommended guide wire	0.035"	Sizes	ø 3.0-12.0 mm; L: 20-250 mm
Tip	Low entry profile, colored	Shaft	5.1-5.4F, dual-lumen, hydrophobic coating
Balloon material	SCP (Semi-Crystalline Polymer), controlled compliance	Usable length	90, 130 and 170 cm
Balloon coating	Hydrophobic patchwork coating	Guide wire lumen	Hydrophobic coating

Balloon ø (mm)				Catheter length 90 cm									
NP ^a RBP ^a				Balloon length (mm)									
				20	40	60	80	100	120	150	170	200	250
5F	3.0	9	21	428777	428786	428795	428804	-	-	-	-	-	-
	4.0	9	20	428778	428787	428796	428805	428814	428823	-	-	-	-
	5.0	9	19	428779	428788	428797	428806	428815	428824	428833	-	428843	428848
	6.0	9	16	428780	428789	428798	428807	428816	428825	428834	-	428844	428849
	7.0	9	15	428781	428790	428799	428808	428817	-	-	-	-	-
6F	8.0	9	14	428782	428791	428800	428809	428818	-	-	-	-	-
	9.0	9	12	428783	428792	428801	428810	428819	-	-	-	-	-
	10.0	8	11	428784	428793	428802	428811	428820	-	-	-	-	-
7F	12.0	8	10	-	428794	428803	428812	428821	-	-	-	-	-

Balloon ø (mm)				Catheter length 130 cm									
NP ^a RBP ^a				Balloon length (mm)									
				20	40	60	80	100	120	150	170	200	250
5F	3.0	9	21	428851	428860	428869	428878	428887	428896	428905	428910	428915	428920
	4.0	9	20	428852	428861	428870	428879	428888	428897	428906	428911	428916	428921
	5.0	9	19	428853	428862	428871	428880	428889	428898	428907	428912	428917	428922
	6.0	9	16	428854	428863	428872	428881	428890	428899	428908	428913	428918	428923
	7.0	9	15	428855	428864	428873	428882	428891	428900	428909	428914	428919	428924
6F	8.0	9	14	428856	428865	428874	428883	428892	428901	-	-	-	-
	9.0	9	12	428857	428866	428875	428884	428893	428902 ^b	-	-	-	-
	10.0	8	11	428858	428867	428876	428885	428894	428903 ^b	-	-	-	-
7F	12.0	8	10	-	428868	428877	428886	428895	428904 ^b	-	-	-	-

				Catheter length 170 cm									
Balloon ø (mm)				NP ^a	RBP ^a	Balloon length (mm)							
				20	40	60	80	100	120	150	170	200	250
5F	4.0	9	20	-	428935 ^b	428944 ^b	428953 ^b	-	428971 ^b	428980 ^b	-	428990 ^b	428995 ^b
	5.0	9	19	428927 ^b	428936 ^b	428945 ^b	428954 ^b	-	428972 ^b	428981 ^b	-	428991 ^b	428996 ^b
	6.0	9	16	428928 ^b	428937 ^b	428946 ^b	428955 ^b	-	428973 ^b	428982 ^b	-	428992 ^b	428997 ^b
	7.0	9	15	428929 ^b	428938 ^b	428947 ^b	428956 ^b	-	-	-	-	-	-
6F	8.0	9	14	428930 ^b	428939 ^b	428948 ^b	428957 ^b	-	-	-	-	-	-
	9.0	9	12	-	428940 ^b	428949 ^b	428958 ^b	-	-	-	-	-	-
	10.0	8	11	-	428941 ^b	428950 ^b	428959 ^b	-	-	-	-	-	-

^aNP = Nominal Pressure and RBP = Rated Burst Pressure in atm (1 atm = 1.013 bar)

^b8 weeks pre-order only

1. BIOTRONIK data on file. Compared to leading competitors 6.0 mm diameter balloon size. IIB(P) 39/217; 2. BIOTRONIK data on file. Compared to leading competitors 6.0x200 mm balloon size. *Note for Australia: Passeo-35 Xeo not approved by Therapeutic Goods Administration for use in the common iliac arteries. **Indication as per IFU. **Passeo and Xeo are trademarks or registered trademarks of the BIOTRONIK Group of Companies.** EverCross is a trademark or registered trademark of the Medtronic Group of companies. Mustang is a trademark or registered trademark of the Boston Scientific Group of companies. Ultraverse is a trademark or registered trademark of Becton, Dickinson and Company or its affiliates. Armada is a trademark or registered trademark of the Abbott Group of companies.

BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronic.com
www.biotronic.com

© 2019 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

BIOTRONIK
excellence for life